


Product Name: Recombinant Human PTC3 Protein
Catalog #: PEH9751



Summary

Name	PTCD3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96EY7
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VTGIEEVVIP KKKTWDKVAV LQALASTVNR DTTAVPYVFQ DDPYLMPASS LESRFLLAK KSGENVAKFI INSYPKYFQK DIAEPHIPCL MPEYFEPQIK DISEAALKER IELRKVKASV DMFDQLLQAG TTVSLETTNS

image



Background

Alternative Names	Pentatricopeptide repeat domain-containing protein 3, mitochondrial (28S ribosomal protein S39, mitochondrial) (MRP-S39) (Mitochondrial small ribosomal subunit protein mS39) (Transformation-related gene 15 protein) (TRG-15)
--------------------------	---

Note

Product Name: Recombinant Human PTCD3 Protein
Catalog #: PEH9751



For Research Use Only , Not for Diagnostic Use