


Product Name: Recombinant Human NGLY1 Protein
Catalog #: PEH9752



Summary

Name	NGLY1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96IV0
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IRDLIAIERS SRLDGSNKSH KVKSSQQPAA STQLPTTPSS NPSGLNQHTR NRQGQSSDPP SASTVAADSA ILEVLSQNIQ HVLVYENPAL QEKALACIPV QELKRKSQEK LSRARKLDKG

image



Background

Alternative Names Peptide-N(4)-(N-acetyl-beta-glucosaminyl)asparagine amidase (PNGase) (hPNGase) (EC 3.5.1.52) (N-glycanase 1) (Peptide:N-glycanase)

Note

For Research Use Only , Not for Diagnostic Use