


Product Name: Recombinant Human FKBP10 Protein
Catalog #: PEH9753



Summary

Name	FKB10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96AY3
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FVRYHYNCSL LDGTQLFTSH DYGAPEATL GANKVIEGLD TGLQGMCVGE RRQLIVPPHL AHGESGARGV PGSVLLFEV ELVSREDGLP TGYLFVWHKD PPANLFEDMD LNKDGEVPPE EFSTFIKAQV SEGKGRMLPG QDPEKTIGDM FQNQDRNQDG KITVDELKLLK SDEDEERVHE

image



Background

Alternative Names	Peptidyl-prolyl cis-trans isomerase FKBP10 (PPIase FKBP10) (EC 5.2.1.8) (65 kDa FK506-binding protein) (65 kDa FKBP) (FKBP-65) (FK506-binding protein 10) (FKBP-10) (Immunophilin FKBP65) (Rotamase)
--------------------------	--

Note

Product Name: Recombinant Human FKBP10 Protein
Catalog #: PEH9753



For Research Use Only , Not for Diagnostic Use