


Product Name: Recombinant Human NUDT12 Protein
Catalog #: PEH9754



Summary

Name	NUDT12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BQG2
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TNEVEECENY FSKTLLDRKS EKRNSDWLL AKESHPATVF ILFSDLNPLV TLGGNKESFQ QPEVRLCQLN YTDIKDYLAQ PEKITLIFLG VELEIKDKLL NYAGEVPREE EDGLVAWFAL GIDPIAAEEF QQRHENCYFL HPPMPALLQL

image



Background

Alternative Names NAD-capped RNA hydrolase NUDT12 (DeNADding enzyme NUDT12) (EC 3.6.1.-) (NADH pyrophosphatase NUDT12) (EC 3.6.1.22) (Nucleoside diphosphate-linked moiety X motif 12) (Nudix motif 12)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human NUDT12 Protein
Catalog #: PEH9754

