


Product Name: Recombinant Human TAMM41 Protein
Catalog #: PEH9758



Summary

Name	TAM41_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96BW9
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PKIITSIQNN YGAGVYYNSL IMCNGRLIKY GVISTNVLIE DLLNWNLYI AGRLQKPVKI ISVNEDVTLR SALDRNLKSA VTA AFLMLPE SFSEEDLFIE IAGLSYSGDF RMVVGEDKTK VLNIVKPNIA HFRELYGSIL QENPQVVYKS QQGWLEIDKS PEGQFTQLMT LPKTLQQQIN HIMDPPGKNR DVEETLFQVA HDPDCGDVVR LGLSAIVRPS SIRQSTKGIF TAGKSFGNPC VTYLLTEWLP HSWLQCKALY LLGACEMLSF DGHKLGYSK VQTGITAAEP GGRTMSDHWQ

image



Background

Alternative Names	Phosphatidate cytidyltransferase, mitochondrial (EC 2.7.7.41) (CDP-diacylglycerol synthase) (CDP-DAG synthase) (Mitochondrial translocator assembly and maintenance protein 41 homolog) (TAM41)
--------------------------	---

Product Name: Recombinant Human TAMM41 Protein
Catalog #: PEH9758



Note

For Research Use Only , Not for Diagnostic Use