


Product Name: Recombinant Human PI4KB Protein
Catalog #: PEH9759



Summary

Name	PI4KB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBF8
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VEVLECFNFD TTSVPAIPE NRIRSTRSVE NLPECGITHE QRAGSFSTVP NYDNDDEAWS VDDIGELQVE LPEVHTNSCD NISQFSVDSI TSQESKEPVF IAAGDIRRRRL SEQLAHTPTA FKRPEDPSA VALKEPWQEK

image



Background

Alternative Names Phosphatidylinositol 4-kinase beta (PI4K-beta) (PI4Kbeta) (PtdIns 4-kinase beta) (EC 2.7.1.67) (NPIK) (PI4K92) (PI4KIII)

Note

For Research Use Only , Not for Diagnostic Use