


Product Name: Recombinant Human PTER Protein
Catalog #: PEH9760



Summary

Name	PTER_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96BW5
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PPCQEAIKKE PIVMKNLYWI QKNAYSHKEN LQLNQETETI KEELLYFKAN GGGALVENTT TGISRDTQTL KRLAETGVH IISGAGFYVD ATHSSETRAM SVEQLTDVLM NEILHGADGT SIKCGIIGEI GCSWPLTESE RQVLQATAHA QAQLGCPVII HPGRSSRAPF QIIRILQEAG ADISKTVMSE LDRITLDKKE LLEFAQLGCI LEYDLFGTEL LHYQLGPDID MPDDNKRIRR VRLVVEEGCE

image



Background

Alternative Names	Phosphotriesterase-related protein (EC 3.1.-.-) (Parathion hydrolase-related protein) (hPHRP)
--------------------------	---

Note

Product Name: Recombinant Human PTER Protein
Catalog #: PEH9760



For Research Use Only , Not for Diagnostic Use