


**Product Name: Recombinant Human PLEKHA1 Protein**  
**Catalog #: PEH9765**



## Summary

<b>Name</b>	PKHA1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9HB21
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	31 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	SQPNSDNLSR HGECGKKQVS YRTDIVGGVP IITPTQKEEV NECGESIDRN NLKRSQSHLP YFTP KPPQDS AVIKAGYCVK QGAVMKNWKR RYFQLDENTI GYFKSELEKE PLRVIPLKEV HKVQECKQSD IMMRDNLFEI VTTSRTFYVQ ADSPEEMHSW IKAVSGAIVA QRGPGRSASS EHPPGPSESK HAFRPTNAAT ATSHSTASRS NSLVSTFTME KRGFYESLAK VKPGNFKVQT VSPREPASKV TEQALLRPQS

## image



## Background

<b>Alternative Names</b>	Pleckstrin homology domain-containing family A member 1 (PH domain-containing family A member 1) (Tandem PH domain-containing protein 1) (TAPP-1)
--------------------------	---

**Product Name: Recombinant Human PLEKHA1 Protein**  
**Catalog #: PEH9765**



---

**Note**

For Research Use Only , Not for Diagnostic Use