


Product Name: Recombinant Human PLEKHD1 Protein
Catalog #: PEH9766



Summary

Name	PLHD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NEE1
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NAQLGEAMIK SLEAQLQLA KEKQEYLDKL MEETEELCLQ REQREELERL NQVLEAEKQQ FEEVVQELRM EQEQIKRELE LTARCLKGVE QEKKELRHLL ESLQQTLEEL SIEKKKTLEM LEENENHLQT LANQSEQPPP SGGLHSNLRQ IEEKMQQLLE EKLLAEKRMK ENEERSALE EEREFYSSQS QALQNSLQEL TAEKQQAERE LKAEVKVRMD LERRLREAEG ALRSLEQGLN SKVRNKEKEE RMRADVSHLK RFFEECIRNA ELEAKMPVIM

image



Background

Alternative Names	Pleckstrin homology domain-containing family D member 1 (PH domain-containing family D member 1)
--------------------------	--

Product Name: Recombinant Human PLEKHD1 Protein
Catalog #: PEH9766



Note

For Research Use Only , Not for Diagnostic Use