



Product Name: Recombinant Human PLEKHH1 Protein
Catalog #: PEH9767



Summary

Name	PKHH1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9ULM0
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DKMQELEQRL LEAEQRAENA ETQVGVMEEK VKLSNLKNVD SEGSLHRKYQ ELLKAIKGD ELISQLEAQL EKQKQMRAEE AKTVQEKA AK IKWVTLKLA KLEMENQHLK SHNQRLVEQV GSLQDALEAI QIAPSRKLLV PPYGAAEQDS VPSEPGIQPM GDDSGSQAQG LKAAVLAPSP GALQSKDSVS EAASPLEDSS SSTVHSGETV

image



Background

Alternative Names Pleckstrin homology domain-containing family H member 1 (PH domain-containing family H member 1)

Note

Product Name: Recombinant Human PLEKHH1 Protein
Catalog #: PEH9767



For Research Use Only , Not for Diagnostic Use