


Product Name: Recombinant Human MTPAP Protein
Catalog #: PEH9769



Summary

Name	PAPD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NVV4
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RSNKQLFELL CYAESIDDQL NTLLKEFQLT EENTKLRYLT CSLIEDMAAA YFPDCIVRPF GSSVNTFGKL GCDLDMFLDL DETRNLSAHK ISGNFLMEFQ VKNVPSERIA TQKILSVLGE CLDHFGPGCV

image



Background

Alternative Names	Poly(A) RNA polymerase, mitochondrial (PAP) (EC 2.7.7.19) (PAP-associated domain-containing protein 1) (Polynucleotide adenylyltransferase) (Terminal uridylyltransferase 1) (TUTase 1) (mtPAP)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MTPAP Protein
Catalog #: PEH9769

