


Product Name: Recombinant Human DHX35 Protein
Catalog #: PEH9772



Summary

Name	DHX35_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H5Z1
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WRPGTEGPGV SISEERQSLA ENSGTTVVYN PYAALSIEQQ RQKLPVFKLR NHILYLIENY QTVVIVGETG CGKSTQIPQY LAEAGWTAEG RVVGVQTQPRR VAAVTVAGRV AEERGAVLGH EVGYCIRFDD CTDQLATRIK FLTDGMLVRE MMVDPLLTKY SVIMLDEAHE RTLYTDIAIG LLKKIQKKRG DLRLIVASAT LDADKFRDFF NQNETSDFPAR

image



Background

Alternative Names Probable ATP-dependent RNA helicase DHX35 (EC 3.6.4.13) (DEAH box protein 35)

Note

Product Name: Recombinant Human DHX35 Protein
Catalog #: PEH9772



For Research Use Only , Not for Diagnostic Use