


Product Name: Recombinant Human MCF2L2 Protein
Catalog #: PEH9774



Summary

Name	MF2L2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86YR7
Host	E.coli
Species	Human
Molecular Mass	58 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VKSEEIFESH HERGNPELEQ QARLGDLSPR RRIIRDLEET EEIYIKEIKS IIDGYITPMD FIWLKHLIPD VLQNNKDFLF GNIRELYEFH NRTFLKELEK CAENPELLAH CFLKRKEDLQ IYFKYHKNLP RARAIWQECQ DCAYFGVCQR QLDHNLPLFK YKGPSQRLLI KYQMLLKGLL DFESPEDMEI DPGELGGSAAK DGPKRTKDSA FSTELQQALA VIEDLIKSCS LAVDLAAVTE CPDDIGKLGK LLLHGPFVSVW TIHKDRYKMK DLIRFKPSQR QIYLFERIV FCKIRMEPGD QGLSPHYFSK KTMKLMTLIS RQLGRGSHRK FEIASRNGLE KYILQAASKE IRDCWFSEIS KLLMEQQNNI KDQGNPQFEM STSKGSGAGS GPWIKNMERA TTSKEDPASS TGGIKGCSSR EFSSMDTFED CEQAEDMEKE SSALSLAGLF QSDDSHETCS SKSAFLERGE SSQGEKEERD EEETATRSTE

image



Background

Product Name: Recombinant Human MCF2L2 Protein
Catalog #: PEH9774



Alternative Names

Probable guanine nucleotide exchange factor MCF2L2 (Dbs-related Rho family guanine nucleotide exchange factor) (MCF2-transforming sequence-like protein 2)

Note

For Research Use Only , Not for Diagnostic Use