


Product Name: Recombinant Human ZC3H12C Protein
Catalog #: PEH9775



Summary

Name	ZC12C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9C0D7
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KPEESQTTSK EAKKPPDVVR EYQTKLEFAL KLGYSSEQVQ LVLNKLGTDA LINDILGELV KLGNKSEADQ TVSTINTITR ETSSLESQRS ESPMQEIVTD DGENLRPIVI DGSNVAMSHG NKEVFSCRGI KLAVDWFLER

image



Background

Alternative Names	Probable ribonuclease ZC3H12C (EC 3.1.-.-) (MCP-induced protein 3) (Zinc finger CCCH domain-containing protein 12C)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use