


Product Name: Recombinant Human PRR35 Protein
Catalog #: PEH9779



Summary

Name	PRR35_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P0CG20
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LQPRGPVPGS PEHVGEDLTR ALGDYARVEQ RLGQLGPAGG LAPRPLREQL GKIRLELLTI HQALEQAVRP PDAPLDLSVK RAPAKGPQAL GEAWGRPELG PVLTTGGTPEP PGMLGPAAPQ PFSGHTTKCE ADSSVPPPGL

image



Background

Alternative Names Proline-rich protein 35 (Uncharacterized protein RJD1)

Note

For Research Use Only , Not for Diagnostic Use