


Product Name: Recombinant Human ABHD8 Protein
Catalog #: PEH9780



Summary

Name	ABHD8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96I13
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VGPLESVSS DGYTFVEVKP GRVLRVKHAG PAPAAAPPPP SSASSDAAQG DLSGLVRCQR RITVYRNGRL LVENLGRAPR ADLLHGQNGS GEPPAALEVE LADPAGSDGR

image



Background

Alternative Names Protein ABHD8 (EC 3.-.-) (Alpha/beta hydrolase domain-containing protein 8) (Abhydrolase domain-containing protein 8)

Note

For Research Use Only , Not for Diagnostic Use