


Product Name: Recombinant Human FAM184A Protein
Catalog #: PEH9787



Summary

Name	F184A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NB25
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RKKLIEDYEG KLNKAQSFYE RELDTLKRSQ LFTAESLQAS KEKEADLRKE FQGQEAILRK TIGKLKTELQ MVQDEAGSLL DKCQKLQTAL AIAENNVQVL QKQLDDAKEG EMALLSKHKE VESELAARE RLQQQASDLV LKASHIGMLQ ATQMTQEVTI KDLESEKSRV NERLSQLEEE RAFLRSKTQS LDEEQKQIL ELEKKVNEAK RTQQEYYERE LKNLQSRLEE EVTQLNEAHS KTLEELAWKH HMAIEAVHSN AIRDKKKLQM DLEEQHNKDK LNLEEDKNQL QQELENLKEV LEDKLNTANQ EIGHLQDMVR KSEQGLGSAE GLIASLQDSQ ERLQNELDLT KDSLKETKDA LLNVEGELEQ

image



Background

Alternative Names	Protein FAM184A
--------------------------	-----------------

Product Name: Recombinant Human FAM184A Protein
Catalog #: PEH9787



Note

For Research Use Only , Not for Diagnostic Use