


## Summary

<b>Name</b>	F184B_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9ULE4
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	38 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	RLQHLSHEA TPQGRLPQES PETKSEPGQG PEMQEVLLEV QRLRVENQQL SKDYARKAEE LQATYERENE AIRQAMQQSV SQALWQWQEK ESDLRKNFQV QESALQAQVR KLEGDLEHRG RKISDLKKYA QKLKERIQDL DVQLKEARQE NSELKGTAKK LGEKLAVAKD RMMLQECRGT QQTDAMKTEL VSENKVLREE NDLEAGNLHP QQDQSCLKEC PCMKGGTDMQ TKKEASAETE YMKQQYEEDL RKIKHQTEEE KKHLKDQLVK RLEDLVKKHT VEIKSVRSSV EAERKKLQRE VEAQLLEVRK

## image



## Background

**Alternative Names** Protein FAM184B

**Product Name: Recombinant Human FAM184B Protein**  
**Catalog #: PEH9788**



---

**Note**

For Research Use Only , Not for Diagnostic Use