



Product Name: Recombinant Human FAM193A Protein
Catalog #: PEH9790



Summary

Name	F193A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P78312
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLQFINSET KPVSSSTRAAK RARHKQRKLE EKARLEAEAR AREHLHLQEE QRRREEEED EEEEDRFKEE FQRLQELQKL RAVKKKKKER PSKDCPKLDM LTRNFQAATE SVPNSGNIHN GSLEQTEEPE TSSHSPSRHM NHSEPRPGLG ADGDAADPVD TRDSKFLLPK EVNGKQHEPL SFFFDIMQHH

image



Background

Alternative Names Protein FAM193A (Protein IT14)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FAM193A Protein
Catalog #: PEH9790

