


Product Name: Recombinant Human FAM47E Protein
Catalog #: PEH9791



Summary

Name	FA47E_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZV65
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HSRQLVFPRK GLDDFRKGCP PCTGLVTQVP VEGFLPQIYH RAPQLAPKKR QIKLLKEADV LSKLSPAQQA RKAFLDVEA HLTPHPLALY LNLEEAMPIE LLSKVLEVLD PDRKLEDTWA YCQDTRKGMK EPTKLLKKHS TQVYLGPSKK TSVSNAGQWL YEEKPHKMDL LHENGPRPGL HENGISDIDE EFILKQFDID YETKPSHDAL HTMKLNQVPL ELKRSVGLSK LQETEFFQKL GYERKLQKPQ NPYKPKWVKM RYGAWYLNPK LWKKQRVDEP LVDPEVSHKA QEENFKKELQ EQEELLADLH GTVAFKDFIL SRGYRTPRFL

image



Background

Alternative Names Protein FAM47E

Product Name: Recombinant Human FAM47E Protein
Catalog #: PEH9791



Note

For Research Use Only , Not for Diagnostic Use