


Product Name: Recombinant Human FAM81A Protein
Catalog #: PEH9795



Summary

Name	FA81A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TBF8
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NKEQQAAKLI LETKIKDAEG QISQLLNRVD LSISEQSTKL KMSHRDSNHQ LQLLDTKFKG TVEELSNQIL SARSWLQREQ ERIEKELLQK IDQLSLIVKE NSGASERDME KKLSQMSARL DKIEEGQKKT FDGQRTRQEE EKMHGTRITKL ELQMNQNIKE MKADEVNAGFT

image



Background

Alternative Names Protein FAM81A

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FAM81A Protein
Catalog #: PEH9795

