


Product Name: Recombinant Human NSMAF Protein
Catalog #: PEH9796



Summary

Name	FAN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92636
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EHDVSVDTIS LNAASTLLVS GTKEGTVNIW DLTTATLMHQ IPCHSGIVCD TAFSPDSRHV LSTGTDGCLN VIDVQTGMLI SSMTSDEPQR CFVWDGNSVL SGSQSGELLV WDLLGAKISE RIQGHTGAVT CIWMNEQCSS IITGGEDRQI

image



Background

Alternative Names	Protein FAN (Factor associated with neutral sphingomyelinase activation) (Factor associated with N-SMase activation)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use