


Product Name: Recombinant Human NKD2 Protein
Catalog #: PEH9806



Summary

Name	NKD2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q969F2
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HSSGSSKTLR VKLTVSPEPS SKRKEGPPAG QDREPTRCRM EGELAEPRV ADRRLSAHVR RPSTDPQPCS ERGPYCVDEN TERRNHLYDL AGIENYTSRF GPGSPPVQAK QEPQGRASHL QARSRSQEPD THAVHHRRSQ VLVEHVVPAS EPAARALDTQ

image



Background

Alternative Names Protein naked cuticle homolog 2 (Naked-2) (hNkd2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human NKD2 Protein
Catalog #: PEH9806

