


Product Name: Recombinant Human NDRG3 Protein
Catalog #: PEH9807



Summary

Name	NDRG3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UGV2
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EIKPLLNDKN GTRNFQDFDC QEHDITTHG VVHVTIRGLP KGNRPVILTY HDIGLNHKSC FNAFFNFEDM QEITQHFAVC HVDAPGQQEG APSFPTGYQY PTMDELAEML PPVLTHLSLK SIIGIGVGAG AYILSRFALN HPELVEGLVL INVDPCAAGW IDWAASKLSG LTTNVVDIIL AHHFGQEELQ ANLDLIQTYR MHIAQDINQD NLQLFLNSYN GRRDLEIERP ILGQNDNKS K TLKCSTLLVV GDNSPAVEAV VECNSRLNPI NTTLLKMADC GGLPQVVQPG

image



Background

Alternative Names	Protein NDRG3 (N-myc downstream-regulated gene 3 protein)
--------------------------	---

Product Name: Recombinant Human NDRG3 Protein
Catalog #: PEH9807



Note

For Research Use Only , Not for Diagnostic Use