


Product Name: Recombinant Human NEDD1 Protein
Catalog #: PEH9808



Summary

Name	NEDD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NHV4
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEENENRDLT AESKKIYMGK QESKDSFKQL AKLVTSGAES GNLNTSPSSN QTRNSEKFEK PENEIEAQLI CEPPINGSST PNPKIASSVT AGVASSLSEK IADSIGNNRQ NAPLTSIQIR FIQNMIQETL DDFREACHRD IVNLQVEMIK QFHMQLNEMH

image



Background

Alternative Names	Protein NEDD1 (Neural precursor cell expressed developmentally down-regulated protein 1) (NEDD-1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human NEDD1 Protein
Catalog #: PEH9808

