


Product Name: Recombinant Human PPP1R21 Protein
Catalog #: PEH9810



Summary

Name	PPR21_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZMI0
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QKLAQEYSKL RAQNQVLKKG VVDEQANSAA LKEQLKMKDQ SLRKLQQEMD SLTFRNLQLA KRVELLQDEL ALSEPRGKKN KKSCESSQL SQEQKSVFDE DLQKKIEENE RLHIQFFEAD EQHKHVEAEL RSRLATLETE AAQHQA VVDG LTRKYMETIE KLQNDKAKLE VKSQTLEKEA KECRLRTEEC QLQLKTLHED LSGRLEESLS IINEKVPFND TKYSQYNALN

image



Background

Alternative Names	Protein phosphatase 1 regulatory subunit 21 (Coiled-coil domain-containing protein 128) (KLRAQ motif-containing protein 1)
--------------------------	--

Note

Product Name: Recombinant Human PPP1R21 Protein
Catalog #: PEH9810



For Research Use Only , Not for Diagnostic Use