


Product Name: Recombinant Human PPP1R27 Protein
Catalog #: PEH9811



Summary

Name	PPR27_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86WC6
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LADRSVRFPN DVLFLDHIRQ GDLEQVGRFI RTRKVSLATI HPSGLAALHE AVLSGNLECV KLLVKYGADI HQRDEAGWTP LHIACSDGYG DIARYLISLG ADRDATNDDG DLPSDLIDPD YKELVELFKG

image



Background

Alternative Names	Protein phosphatase 1 regulatory subunit 27 (Dysferlin-interacting protein 1) (Toonin)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use