


Product Name: Recombinant Human SEC23A Protein
Catalog #: PEH9817



Summary

Name	SC23A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q15436
Host	E.coli
Species	Human
Molecular Mass	91 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ISELNQPAEL LPQFSSIEYV VLRGPQMPLI FLYVVDTCTME DEDLQALKES MQMSLSLLPP TALVGLITFG RMVQVHELGC EGISKSIVFR GTKDLSAKQL QEMLGSLKVP LTQATRGPQV QQPPPSNRFL QPVQKIDMNL TDLLGELQRD PWPVPQGKRP LRSSGVALSI AVGLLECTFP NTGARIMMFI GGPATQGPGM VVGDELKTPI RSWHDIDKDN AKYVKKGTKH FEALANRAAT TGHVIDIYAC ALDQTGLLEM KCCPNLTGGY MVMGDSFNFS LFKQTFQRVF TKDMHGQFKM GFGGTLEIKT SREIKISGAI GPCVSLNSKG PCVSENEIGT GGTCQWKICG LSPTTTLAIY FEVVNQHNAP IPQGGRGAIQ FVTQYQHSSG QRRIRVTTIA RNWADAQTQI QNIAASFDQE AAAILMARLA IYRAETEEGP

image



Background

Product Name: Recombinant Human SEC23A Protein
Catalog #: PEH9817



Alternative Names

Protein transport protein Sec23A (hSec23A) (SEC23-related protein A)

Note

For Research Use Only , Not for Diagnostic Use