


Product Name: Recombinant Human ZBED8 Protein
Catalog #: PEH9820



Summary

Name	ZBED8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IZ13
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LEIAQIVLGP DAQKKLQQVP LSDDVIHSRI DEMSQDILQQ VLEDIKASPL KVGIQLAETT DMDDCSQLMA FVRYIKEREI

image



Background

Alternative Names	Protein ZBED8 (Transposon-derived Buster3 transposase-like protein) (Zinc finger BED domain-containing protein 8)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use