


Product Name: Recombinant Human NBPF5P Protein
Catalog #: PEH9825



Summary

Name	NBPF5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86XG9
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ERAEMNILEI NQELRSQLAE SNQQFRDLKE KFLITQATAY SLANQLKKYK CEEYKDIIDS VLRDELQSME KLAELRQAE ELRQYKALVH SQAQELTQLR EKLREGRDAS RWLNKHLKTL LTPDDPKSQ GQDLREQLAE GHRLAEHLVH KLSPENDEDE DEDEDDKDEE VEKVQESPAP REVQKTEEKE VPQDSLEECA VTCSNSHNPS NSNQPHRSTK ITFKEHEVDS ALVVESEHPH DEEEEALNIP PENQNDHEEE EGKAPVPPRH HDKSNSYRHR EVSFLALDEQ KVCSAQDVAR DYSNPKWDET SLGFLDTPLA RRESVALKGR

image



Background

Alternative Names	Putative neuroblastoma breakpoint family member 5
--------------------------	---

Product Name: Recombinant Human NBPF5P Protein
Catalog #: PEH9825



Note

For Research Use Only , Not for Diagnostic Use