


Product Name: Recombinant Human PRORS1P Protein
Catalog #: PEH9827



Summary

Name	PRXD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NEY8
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAGAEALGAAL EQRLGALAIH TEVVEHPEVF TVEEMMPHIQ HLKGAHSKNL FLKDKKKKNY WLVTVLHDRQ INLNELAKQL GVGSGNLRFA DETAMLEKLK VGQGCATPLA LFCDGGDVKF VLDSAFLEGG HEKVYFHPMT NAATMGLSPE DFLTFVKMTG HDPIILNFD

image



Background

Alternative Names	Putative prolyl-tRNA synthetase associated domain-containing protein 1 (PrdX deacylase domain-containing protein 1) (Prolyl-tRNA synthetase associated domain-containing protein 1 pseudogene)
--------------------------	--

Note



Product Name: Recombinant Human PRORS1P Protein
Catalog #: PEH9827

For Research Use Only , Not for Diagnostic Use