


Product Name: Recombinant Human MAPKAPK5-AS1 Protein
Catalog #: PEH9831



Summary

Name	MAAS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N8E1
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ILSGAPTAGR GGGCSAALSP RGRGSKGLGT RAPGPRGDGQ QPPLGTGGDE DPGAGSASAG GSRLAAAAAA EAAAPGDRSL CWAPRGRLA SSPAGEAGSW GRARRGGPRP GAPCKGLAGP

image



Background

Alternative Names	Putative uncharacterized protein encoded by MAPKAPK5-AS1 (MAPKAPK5 antisense RNA 1) (MAPKAPK5 antisense gene protein 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use