


Product Name: Recombinant Human RABGAP1L Protein
Catalog #: PEH9835



Summary

Name	RBG10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	B7ZAP0
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FEERENRRLQ EASMRLEQEN DDLAHELVT S KIALRNDLDQ AEDKADVLNK ELLTKQRLV ETEEEKRKQE EETAQLKEVF RKQLEKAEYE IKKTTAIIAE YKQICSQLST RLEKQQAASK EELEVVKGKM MACKHCSDIF SKEGALKLAA TGREDQGIET DDEKDSLKKQ LREMELELAQ TKLQLVEAKC KIQELEHQRG ALMNEIQAAK NSWFSKTLNS

image



Background

Alternative Names Rab GTPase-activating protein 1-like, isoform 10

Note

Product Name: Recombinant Human RABGAP1L Protein
Catalog #: PEH9835



For Research Use Only , Not for Diagnostic Use