


Product Name: Recombinant Human RANBP3 Protein
Catalog #: PEH9840



Summary

Name	RANB3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H6Z4
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSDANRENAA AESGSESSQ EATPEKESLA ESAAAYTKAT ARKCLLEKVE VITGEEAESN VLQMCKLFV FDKTSQSWVE RGRGLRLND MASTDDGTLQ SRLVMRTQGS LRLILNTKLW AQMQIDKASE KSIRITAMDT EDQGVKVFLI SASSKDTGQL YAALHHRILA LRSRVEQEQE AKMPAPEPGA APSNEEDDSD DDDLAPSGA

image



Background

Alternative Names	Ran-binding protein 3 (RanBP3)
--------------------------	--------------------------------

Note

Product Name: Recombinant Human RANBP3 Protein
Catalog #: PEH9840



For Research Use Only , Not for Diagnostic Use