


Product Name: Recombinant Human RRAGC Protein
Catalog #: PEH9842



Summary

Name	RRAGC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HB90
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MEALTRLHIT VSKAYKVNPD MNFEVFIHKV DGLSDDHKIE TQRDIHQAN DDLADAGLEK LHLSFYLTISI YDHSIFEAFS KVVQKLIPQL PTLENLLNIF ISNSGIEKAF LFDVVSKIYI ATDSSPVDMQ SYELCCDMID VVIDVSCIYG

image



Background

Alternative Names Ras-related GTP-binding protein C (Rag C) (RagC) (EC 3.6.5.-) (GTPase-interacting protein 2) (TIB929)

Note

For Research Use Only , Not for Diagnostic Use