


Product Name: Recombinant Human RIMS3 Protein
Catalog #: PEH9847



Summary

Name	RIMS3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UJD0
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IAVEMRSRVT RQGSRESTDG STNSNSSDGT FIFPTTRLGA ESQFSDFLDG LGPAQIVGRQ TLATPPMGDV HIAIMDRSGQ LEVEVIEARG LTPKPGSKSL PATYIKVYLL ENGACLAKKK TKMTKKTCDP LYQQALLFDE GPQGKVLQVI VWGDYGRMDH KCFMGMAQIM LDELDLSAAV TGWYKLFPTS SVADSTLGSL

image



Background

Alternative Names	Regulating synaptic membrane exocytosis protein 3 (Nim3) (RIM3 gamma) (Rab-3-interacting molecule 3) (RIM 3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RIMS3 Protein
Catalog #: PEH9847

