


Product Name: Recombinant Human ARHGEF37 Protein
Catalog #: PEH9854



Summary

Name	ARH37_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A1IGU5
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDALGRTSNQ LRSFQETFEK VQPPPTTQPL LPGSERQVQA LLSRYGPGKL YQVTSNISGT GTLDTLPRG QIVAILQNKD TKGNSGRWLV DTGGHRGYVP AGKLQLYHVV PSAEELRRQA GLNKDPRCLT

image



Background

Alternative Names Rho guanine nucleotide exchange factor 37

Note

For Research Use Only , Not for Diagnostic Use