


Product Name: Recombinant Human ARHGEF38 Protein
Catalog #: PEH9855



Summary

Name	ARH38_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NXL2
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DIHRSKLLST YSAEELYQAK RKCNDATQEYD INLLEGDLVA VIEQKDPLGS TSRWLVDTGN VKGYVYSSFL KPYNPAKMQK VDAENRFCDD DFENISLFVS SRPASDSVTG TSESSIGDSS SLSGTCGKF ETNGTDVDSF

image



Background

Alternative Names Rho guanine nucleotide exchange factor 38

Note

For Research Use Only , Not for Diagnostic Use