


Product Name: Recombinant Human RNF141 Protein
Catalog #: PEH9860



Summary

Name	RN141_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WVD5
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GSLTYEEFLG RVAELNDVTA KVASGQEKHL LFEVQPGSDS SAFWKVVVRV VCTKINKSSG IVEASRIMNL YQFIQLYKDI TSQAAGVLAQ SSTSEEPDEN SSSVTSCQAS LWMGRVKQLT DEEECCICMD GRADLILPCA

image



Background

Alternative Names RING finger protein 141 (Zinc finger protein 230)

Note

For Research Use Only , Not for Diagnostic Use