



Product Name: Recombinant Human RBM12B Protein
Catalog #: PEH9863



Summary

Name	RB12B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IXT5
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QEEGVNFSVM SSEKMQARSQ SRERGDHSHL FDSKDPPIYS VGAFENFRHQ LEDLRQLDNF KHPQRDFRQP DRHPPEDFRH SSEDFRFPPE DFRHSPEDFR RPREEDFRRP SEEDFRRPWE EDFRRPPEDD FRHPREEDWR RPLEEDWRRP LEEDFRRSPT EDFRQLPEED FRQPPEEDLR WLPEEDFRRP PEEDWRRPPE

image



Background

Alternative Names RNA-binding protein 12B (RNA-binding motif protein 12B)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RBM12B Protein
Catalog #: PEH9863

