


Product Name: Recombinant Human RBM14 Protein
Catalog #: PEH9864



Summary

Name	RBM14_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96PK6
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AGQPSAAYLS MSQGAVANAN STPPPYERTR LSPPRASYDD PYKKAVAMSK RYGSDRRLAE LSDYRRLSES QLSFRRSPTK

image



Background

Alternative Names RNA-binding protein 14 (Paraspeckle protein 2) (PSP2) (RNA-binding motif protein 14) (RRM-containing coactivator activator/modulator) (Synaptotagmin-interacting protein) (SYT-interacting protein)

Note

For Research Use Only , Not for Diagnostic Use