


Product Name: Recombinant Human RALYL Protein
Catalog #: PEH9866



Summary

Name	RALYL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86SE5
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DYHGRVPPPP RAVIPLKRPR VAVTTTRRGK GVFSMKGGSR STASGSTGSK LKSDDELQTIK KELTQIKTKI DSSLGRLEKI EKQQKAEAEA QKKQLEESLV LIQEECVSEI ADHSTEEP AE GGP DADGEEM TDGIEEDFDE

image



Background

Alternative Names RNA-binding Raly-like protein (hRALYL) (Heterogeneous nuclear ribonucleoprotein C-like 3) (hnRNP core protein C-like 3)

Note

For Research Use Only , Not for Diagnostic Use