


Product Name: Recombinant Human SEPHS2 Protein
Catalog #: PEH9870



Summary

Name	SPS2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99611
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QALGLSPSWR LTGFSGMKGU GCKVPQEALL KLLAGLTRPD VRPPLGRGLV GGQEEASQEA GLPAGAGPSP TFPALGIGMD SCVIPLRHGG LSLVQTTDF YPLVEDPYMM GRIACANVLS DLYAMGITEC DNMLMILLSVS QSMSEEEREK VTPLMVKGFR DAAEEGGTAV TGGQTVVNPW IIIGGVATV CQPNEFIMPD SAVVGDLVVL TKPLGTQVAV NAHQWLDNPE RWNKVKMVVS

image



Background

Alternative Names Selenide, water dikinase 2 (EC 2.7.9.3) (Selenium donor protein 2) (Selenophosphate synthase 2)

Note

Product Name: Recombinant Human SEPHS2 Protein
Catalog #: PEH9870



For Research Use Only , Not for Diagnostic Use