


Product Name: Recombinant Human SNTN Protein
Catalog #: PEH9871



Summary

Name	SNTAN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NMZ2
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MGGCMHSTQD KSLHLEGDPN PSAAPTSTCA PRKMPKRISI SKQLASVKAL RKCSBLEKAI ATTALIFRNS SDDSGKLEKA IAKDLLQTQF RNFAEGQETK PKYREILSEL DEHTENKLDLDF EDFMILLLSI TVMSDLLQNI RNVKIMK

image



Background

Alternative Names Sentan (Protein S100-A1-like) (S100 calcium-binding protein A1-like)

Note

For Research Use Only , Not for Diagnostic Use