


**Product Name: Recombinant Human SRR Protein**  
**Catalog #: PEH9874**



## Summary

<b>Name</b>	SRR_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9GZT4
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	37 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	DVEKAHINIR DSIHLTPVLT SSILNQLTGR NLFFKCELFQ KTGSKIRGA LNAVRSLVPD ALERKPKAVV THSSGNHGQA LTYAAKLEGI PAYIVVPQTA PDCKKLAIQY YGASIVYCEP SDESRENVAK RVTEETEGIM VHPNQEPAVI AGQGTIALEV LNQVPLVDAL VVPVGGGGML AGIAITVKAL KPSVKVYAAE PSNADDCYQS KLK GKLM PNL YPPETIADGV KSSIGLNTWP IIRDLVDDIF TVTEDEIKCA TQLVWERMKL LIEPTAGGVV AAVLSQHFQT VSPEVKNICI VLSGGNVDLT SSITWVKQAE RPASYQSVSV

## image



## Background

**Alternative Names** Serine racemase (EC 5.1.1.18) (D-serine ammonia-lyase) (D-serine dehydratase) (EC 4.3.1.18) (L-serine ammonia-lyase) (L-serine dehydratase)

**Product Name: Recombinant Human SRR Protein**  
**Catalog #: PEH9874**



---

(EC 4.3.1.17)

**Note**

For Research Use Only , Not for Diagnostic Use