


Product Name: Recombinant Human SGK3 Protein
Catalog #: PEH9879



Summary

Name	SGK3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96BR1
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ESCPSVSIPS SDEHREKKR FTVYKVLVSV GRSEWVFVRR YAEFDKLYNT LKKQFPAMAL KIPAKRIFGD NFDPDFIKQR RAGLNEFIQN LVRYPELYNH PDVRAFLQMD SPKHQSDPSE DEDERSQKL HSTSQNINLG

image



Background

Alternative Names	Serine/threonine-protein kinase Sgk3 (EC 2.7.11.1) (Cytokine-independent survival kinase) (Serum/glucocorticoid-regulated kinase 3) (Serum/glucocorticoid-regulated kinase-like)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SGK3 Protein
Catalog #: PEH9879

