


**Product Name: Recombinant Human VRK1 Protein**  
**Catalog #: PEH9880**



## Summary

<b>Name</b>	VRK1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q99986
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	46 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MPRVKAAQAG RQSSAKRHLA EQFAVGEIIT DMAKKEWKVG LPIGQGGFGC IYLADMNSSE SVGSDAPCVV KVEPSDNGPL FTELKFYQRA AKPEQIQKWI RTRKLYLGV PKYWGSGLD KNGKSYRFMI MDRFGSDLQK IYEANAKRFS RKTVLQSLR ILDILEYIHE HEYVHGDIKA SNLLNYKNP DQVYLVDYGL AYRYCPEGVH KEYKEDPKRC HDGTIEFTSI DAHNGVAPSR RGDLEILGYC MIQWLTGHLP WEDNLKDPKY VRDSKIRYRE NIASLMDKCF PEKNKPGEIA KYMETVKLLD YTEKPLYENL RDILLQGLKA IGSKDDGKLD LSVVENGGLK AKTITKRRKK EIEESKEPGV EDTEWSNTQT EEAIQTRSRT RKRIVQK

## image



## Background

**Alternative Names** Serine/threonine-protein kinase VRK1 (EC 2.7.11.1) (Vaccinia-related kinase)

**Product Name: Recombinant Human VRK1 Protein**  
**Catalog #: PEH9880**

---



1)

**Note**

For Research Use Only , Not for Diagnostic Use