


Product Name: Recombinant Human ANKRD52 Protein
Catalog #: PEH9883



Summary

Name	ANR52_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NB46
Host	E.coli
Species	Human
Molecular Mass	55 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FCLEFLDNG ADPSLRDRQG YTAVHYAAAY GNRQNLELLL EMSFNCLEDV ESTIPVSPLH LAAYNGHCEA LKTLAETLVN LDVRDHKGRT ALFLATERGS TECVEVLTAA GASALIKERK RKWTPLHAAA ASGHTDSLHL LIDSGERADI TDVMDAYGQT PLMLAIMNGH VDCVHLLLEK GSTADAADLR GRTALHRGAV TGCEDCLAAL LDHDAFVLCR DFKGRTPHIL ASACGHTAVL RTLLQAALST DPLDAGVDYS GYSPMHWASY TGHEDCLELL LEHSPFSYLE GNPFTPLHCA VINNQDSTTE MLLGALGAKI VNSRDAKGRT PLHAAAFADN VSGLRMMLLQH QAEVNATDHT GRTALMTAAE NGQTAAVEFL LYRGKADLTV LDENKNTALH LACSKGHEKC ALMILAETQD LGLINATNSA LQMPHIAAR NGLASVVQAL LSHGATVLAV DEEGHTPALA CAPNKDVADC LALILSTMKP

image



Background

Product Name: Recombinant Human ANKRD52 Protein
Catalog #: PEH9883



Alternative Names

Serine/threonine-protein phosphatase 6 regulatory ankyrin repeat subunit C (PP6-ARS-C) (Serine/threonine-protein phosphatase 6 regulatory subunit ARS-C) (Ankyrin repeat domain-containing protein 52)

Note

For Research Use Only , Not for Diagnostic Use