


Product Name: Recombinant Human SERPINB13 Protein
Catalog #: PEH9884



Summary

Name	SPB13_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9UIV8
Host	E.coli
Species	Human
Molecular Mass	71 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEKEVIENTE AVHQQFQKFL TEISKLTNDY ELNITNRLFG EKTYLFLQKY LDYVEKYYHA SLEPVDFVNA ADESRKKINS WVESKTNEKI KDLFPDGSIS SSTKLVLVNM VYFKGQWDRE FKKENTKEEK FWMNKSTSXS VQMMTQSHSF SFTFLEDLQA KILGIPYKNN DLSEFVLLPN DIDGLEKIID KISPEKLVEW TSPGHMEERK VNLHLPRFEV EDGYDLEAVL AAMGMGDAFS

image



Background

Alternative Names Serpin B13 (HaCaT UV-repressible serpin) (Hurpin) (Headpin) (Peptidase inhibitor 13) (PI-13) (Proteinase inhibitor 13)

Note

Product Name: Recombinant Human SERPINB13 Protein
Catalog #: PEH9884



For Research Use Only , Not for Diagnostic Use