


Product Name: Recombinant Human SERPINB13 Protein
Catalog #: PEH9885



Summary

Name	SPB13_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9UIV8
Host	E.coli
Species	Human
Molecular Mass	84 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LGFDLFKELK KTNDGNIFFS PVGILTAIGM VLLGTRGATA SQLEEVFHSE KETKSSRIKA EEKEVIENTE AVHQQFQKFL TEISKLNDY ELNITNRLFG EKTYLFLQKY LDYVEKYYHA SLEPVDVNA ADESRKKINS WVESKTNEKI KDLFPDGSIS SSTKLVLVNM VYFKGQWDRE FKKENTKEEK FWMNKSTS VQMMTQSHSF SFTFLEDLQA KILGIPYKNN DLSEMFVLLPN DIDGLEKI KISPEKLVW TSPGHMEERK VNLHLPRFEV EDGYDLEAVL AAMGMGDAFS EHKADYSGMS SGSGLYAQKF LHSSFVAVTE EGTEAAAATG IGFTVTSAPG

image



Background

Alternative Names Serpin B13 (HaCaT UV-repressible serpin) (Hurpin) (Headpin) (Peptidase inhibitor 13) (PI-13) (Proteinase inhibitor 13)

Product Name: Recombinant Human SERPINB13 Protein
Catalog #: PEH9885



Note

For Research Use Only , Not for Diagnostic Use